You are not alone!

Join Us!

A community of women
for support and advice.

Dr. Amanda Olson, DPT, PRPC

The Essentials

Perfect with Kegels

VELVET ROSE
Personal Lubricant

ENCHANTED ROSE
Organic Feminine Balm

Shop @:
IntimateRose.com

Manufactured For: Intimate Rose
North Kansas City, MO 64116
Made in PRC

Kegel Exercise System

Includes French, German, Italian, Polish,
Portuguese, Spanish, Swedish, Dutch, and
Turkish Translations.


Keys to Success

The keys to better pelvic floor strength are consistency and progressive resistance. Without consistently training your pelvic floor muscles, you will lose the strength and control you've gained. Without using progressively heavier weights, it is difficult to make significant improvements. The goal is to spend 10-15 minutes every day using one of the weights while being active. While you are up and moving, your pelvic floor muscles will be actively contracting to keep the weight in place. Many users find it convenient to add these weights to their daily morning routine when they wake up, shower, and get ready for the day.

Your progress will depend on multiple factors including: initial pelvic floor strength, age, motivation, adherence to the program, overall pelvic health and history.

Indications for Use

The Intimate Rose Kegel Exercise System is indicated for the strengthening of the perineal pelvic floor muscles by providing resistance to an individual's voluntary contractions of these muscles. It seeks to correct, through exercise, low tone in the pelvic floor, which can cause or contribute to health issues including urinary incontinence and sexual dysfunction. Each color is a different weight. The weights get heavier as the shades get darker.


How To Start

The first step is to find your starter weight. Wash the weights with warm water and soap. Begin with the ivory colored weight. Position yourself in a way that is comfortable with one leg raised on a step or chair. Gently insert the ivory weight by grasping the end of the weight that has the cord attached. Use a water-based lubricant as needed. The cord should remain outside of your body. As you tighten your pelvic muscles around the weight, you may notice the weight naturally pulls up away from your fingers as it is drawn into your vagina. Hold the weight in for two minutes using only your pelvic muscles.

If you successfully hold it in for two minutes, remove the weight and repeat this process with the next darker shade. Continue this process until you have found a challenging, yet comfortable, weight. When you find a weight you cannot hold, use the previous lighter shade as your starter weight. If you have trouble holding the ivory weight, check IntimateRose.com/how for information on building a base level of strength.

Once you have your starter weight identified, your goal is to use this weight for 10-15 minutes one or two times per day. Once you are able to hold the weight for 15 minutes while walking, climbing stairs, or coughing, you are ready for the next darker shade.

Cleaning

It is important to clean the weights before and after each use. Use a mild soap and warm water. Rinse the weights thoroughly. Examine the device for any residual soiling. If found, repeat cleaning and rinsing to remove. Examine for any cracks or additional residual soil. If found, discontinue use and dispose of the weight in the trash. Contact us at Support@IntimateRose.com for a replacement.

Strength Maintenance

If you've reached a sufficient level of pelvic floor strength, it's important to continue a maintenance exercise program. In order to maintain a high level of strength, continue to use the heaviest weight you can for 15 minutes per session three to five times per week.

1. Clean your weights before and after each use.
2. Use the weights up to 15 minutes each session and at most twice per day.
3. Empty your bowel and bladder before each session.
4. If you use a lubricant, be sure it's water based.
5. Consistency is key. Those who add the exercises to their daily routine have the best results!
5. If you're having trouble keeping the weight in, be sure you're not flexing your abdominal muscles too much. Flexing your abs can push the weight out.
6. If you're having trouble finding your pelvic floor muscles, you can attempt to stop the flow of urine when urinating. The feeling of contracting those muscles to stop the urine flow is the same feeling you should have when trying to keep the weight in.

Warnings & Disclaimer

The information included with Intimate Rose products is not a substitute for medical care or advice. If you have, or suspect you have, a health problem, you should always consult a health care professional. You should check with your health care professional and completely read all documentation before using any medical device.

- Do not use during the first six weeks post-partum, after pelvic

surgery, or during pregnancy without clearance from a health care professional.

- Do not use with grade 4 pelvic organ prolapse.
- Do not use with a vaginal or urinary infection.
- Remove other intra-vaginal devices such as tampons, pessaries, or e-strings before using this product.
- Do not use if you experience pain while inserting the device.
- Do not share your weights with others.
- Do not use if there are any visible cracks in the silicone.

Carefully following use and care instructions will extend the life of your devices. Before each use, check for signs of wear or damage including splits/cracks/pitting in silicone, cracks/damage to tail, etc. Pull down on the tail and examine it carefully for signs of potential damage to prevent any difficulty removing the device. If signs of wear or damage are found, discontinue use and dispose in the trash.

Contact **Support@IntimateRose.com** for a replacement.

The weight should be easy to remove with little or no resistance. To remove, relax your pelvic floor muscles and pull gently on the silicone string. If you encounter any resistance removing the device, apply a generous amount of water based lubricant, perform a squat, or raise one leg on a step to help aid removal.

What to Expect

Expect your pelvic muscles to be sore and tired after using the weights, especially in the beginning stages. If extreme soreness occurs, take a day or two off. As you continue to get stronger, the amount of discomfort will decrease. If you have serious vaginal pain, irritation, inflammation, or infection, discontinue your kegel exercise routine.

With regular cleaning and following the guidelines in this manual, you should not encounter any adverse effects from using these weights. If you have questions or concerns, feel free to contact us at **Support@IntimateRose.com**.

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